



## ORIGINAL RESEARCH ARTICLE

## A Mathematical Method to Easily and Accurately Measure the Tumor Volume of Prostate Cancer on Pathological Slides

John Zihou Chen<sup>1</sup>, Monika Wrzolek<sup>2</sup>, Lynne Opitz<sup>2</sup> and Jianli Chen<sup>2\*</sup>

<sup>1</sup>Jericho High Senior School, Jericho, USA

<sup>2</sup>Department of Pathology and Laboratory Medicine, Staten Island University Hospital, USA

\*Corresponding author: Jianli Chen, MD, PhD, Department of Pathology and Laboratory Medicine, Staten Island University Hospital, Northwell Health, NY 10305, USA, Tel: 718-226-4668, Fax: 718-226-1957



### Abstract

**Objective:** For decades, it has been a big challenge to directly and accurately measure an irregularly shaped area for a tumor volume evaluation of prostate cancer. All current methods are not ideal for different reasons. An objective, practical, accurate, and reliable method for measuring of prostate cancer volume is needed, in particular in modern precision medicine setting, pursuing accuracy and quality of health care.

**Methods and materials:** In this study, we applied a mathematical method - the Gauss area calculation formula to the evaluation of irregular shape of prostate cancer areas on H&E slides. To simplify the calculations, we used JavaScript program to achieve that practical easiness.

**Results:** By using this mathematical method, we were able to easily obtain the accurate tumor volume percentage, which has never been so objective, accurate, and practical before. The method can be easily applied using scanned H&E pathological slides - including surgical slides, fine needle aspirate (FNA) slides, and even digitally scanned pictures, particularly when there are separate areas of interest on one slide.

**Conclusion:** To our knowledge, this is the first reported practical solution of this methodological problem. In this project, with a mathematical method as the principle, we developed a JavaScript program to simplify the calculation so the pathologists can easily apply it to obtain an objective and accurate information for the tumor volume of prostate cancer.

### Keywords

Prostate cancer, Tumor volume, Mathematical method, Gauss's area calculation formula

### Introduction

Prostate cancer is the second leading cause of cancer death in American men, behind lung cancer. The American Cancer Society's estimates about 174,650 new diagnoses and about 31,620 deaths from prostate cancer in the United States for 2019 [1]. There are two most significant factors which affect patients' prognosis: "quantity" and "quality". The "quantity" means the tumor volume. Tumor volume shares something in common with staging but with some differences. By the current College of American Pathologists (CAP) prostate cancer protocol, stage 2 prostate cancer is defined as organ confined. This is a very broad category, with tumor volume that may range from 0.1% through 99.9%, as long as the tumor is still confined within the prostate gland. Under this criteria, patients at the same stage but with different tumor volume may have different prognosis [2]. This is particularly true if the tumor volume is large enough to extend to other tissue or organs. On the other hand, even with the same Gleason score, grade group, histological type, and staging, the management of patient with only 1% of tumor volume should not be the same as the one with 99% of tumor volume. Secondly, the "quality", meaning the tumor Gleason score, grade group, and even histology type. All these factors are important information for patient management [3,4] and prognosis [5,6]. This is why the pathologists are required to report the Gleason score, grade group, and percentage of higher pattern (4 or 5) in detail [2] (see Table 1) in the pathological report.

**Table 1:** The Current CAP protocol for prostate cancer grading and Gleason score standard.

Grade group	Gleason score	Definition
1	Less than or equal to 6	Only individual discrete well-formed glands
2	3 + 4 = 7	Predominantly well-formed glands with lesser component of poorly formed/fused/cribriform glands
3	4 + 3 = 7	Predominantly poorly formed/fused/cribriform glands with lesser component (#) of well-formed glands
4	4 + 4 = 8	Only poorly formed/fused/cribriform glands
	3 + 5 = 8	Predominantly well-formed glands and lesser component (##) lacking glands (or with necrosis)
	5 + 3 = 8	Predominantly lacking glands (or with necrosis) and lesser component (##) of well-formed glands
5	9 - 10	Lack gland formation (or with necrosis) with or without poorly formed/fused/cribriform glands (#)

The most recent CAP explanatory notes (posted in August 2019) mentioned that studies have shown prostate cancer volume is a prognostic factor, although the current data are conflicting as to its independent prognostic significance. There are many methods of estimating the amount of tumor in pathology prostate specimens [7-9]. Some reports even mentioned tumor volume can be “eyeballed” by simple visual inspection [10] or it “may be possible” to measure with at least 2 dimensions [11]. Obviously all these methods are highly inaccurate. Perera, et al. [12] compared the contemporary methods and concluded that a reliable good method is still needed.

According to the CAP statement, our basic question is if we were lack of an objective standard method showing accurate tumor volume, how could we compare the data with each other to reach the conclusion that data are conflicting or not? Based on common sense, it is almost absolutely true that prostate cancer volume is a prognostic factor if we simply go to the extremes by comparing 99% with 1% of tumor volume, even though by current CAP staging criteria, both volumes are stage 2, no matter what the Gleason score, grade group, or histology type is.

In order to find out if the differences are significant, for example between 1% versus 5% of tumor volume, we definitely need an accurate tumor volume measurement. We present here a method which allows to easily and accurately measure the tumor volume percentage from the microscopic slides, using a mathematical principle and computerization of the complicated calculations into a JavaScript program. To our knowledge, this is the first reported method which can objectively and directly measure the real tumor volume. We believe this long-term methodological problem can now be marked with a “The end”. The aspects for the potential clinical significance are also discussed.

## Material and Method

The Gauss area formula, a mathematical method to calculate the area of an n-sided, non-self-intersecting

polygon given the coordinates of its vertices, was employed to calculate the tumor area and tissue area separately and then to determine the percentage of tumor area(s) over tissue area.

To approximate areas of tumor and tissue, ordered points were plotted (counterclockwise or clockwise) along the perimeter of the tissue and tumor on an uploaded computer image (JPG, PNG, etc), and the area of resulting irregular polygons found by connecting them was calculated with Gauss’ formula in a JavaScript program: For ordered points  $(x_1, y_1), (x_2, y_2), (x_3, y_3), \dots, (x_n, y_n)$ . The mathematical formula is:

$$Area = \frac{1}{2} \left| \sum_{i=1}^{n-1} x_i y_{i+1} - \sum_{i=1}^{n-1} x_{i+1} y_i \right|$$

Image size and changing dimensions by stretching will not affect tumor area to tissue ratio. JavaScript program steps are as follows: see the [Supplemental Data](#) attached.


## Results

Based on the basic mathematical Gauss’s area calculation formula as principle, we performed the underlying complicated calculation in a JavaScript program. We then built a simple software incorporating the whole process, to obtain the final tumor percentage. The software is available at the following website:

JavaScript program html:

<https://area-of-shape-calculation.johnchen7.repl.co>

To use it, first open the linkage in a window or mac page. And then follow the steps below:



1. Upload a histological H&E picture or any scanned slide;
2. Click the "Tissue" button and start dotting along the perimeter of the tissue; at the last dotting, click the "EndPath" button to close the area which you choose as tissue area;
3. Click the "Tumor" button and then dot around the perimeter of the tumor area; at the last dotting, click the "EndPath" button to close the area which you choose as tumor area;
4. Click the "Result" button and then the percentage of tumor over tissue will be automatically reported out!


Now let's use a few examples to better understand how to easily use it:

#### Example 1:

1. Load the slide by clicking the "Choose File" button.




2. Click the "Tissue" button and then dot the tissue area along the perimeter of the tissue on slide; it does not matter wherever you start. Click the "EndPath" button to end the circling of the tissue area.




3. Click the "Tumor" button and then dot the tumor area along the perimeter of tumor area; click the "EndPath" button to end the circling of the tumor area.



4. Click the "Result" button to get the accurate percentage of tumor area: 23.49%!




#### Example 2: Fine needle aspirate (FNA) slide: 12.99%!



For this FNA slide, our result is 12.99%, which reflects not only the length of the tumor area, but also the actual width. In our daily practice, we usually only measure the length, by which strictly speaking it is not accurate since the tumor may be only partial width of the length. Hence our method is more accurate and the result objectively stands for the real tumor area.

Example 3: If on one slide or picture there are more than one tumor areas, do not worry. This method still works well. Just click "Tumor" button again after finishing the first tumor area. This example includes two irregular tumor areas, with the total percentage 28.13%!



## Discussion

The determination of the tumor percentage based on distribution of the irregularly shaped prostate cancer area(s) in the microscopic slides has a potential clinical significance, with impact on management and prognosis. However, an accurate and practical method of tumor volume measurement has always been challenging, in particular in the prostate cancer which may often be multifocal, with irregularly shaped cancer foci. In the presented project, with a mathematical method as the principle, we developed a simple computerized JavaScript program to simplify the calculation to obtain an objective and accurate percentage of the tumor volume. To our knowledge this is the first paper to reach this practical level and solve this methodological problem. The more accurate results of tumor volume measurement may hopefully allow for reproducible stratification and sub-staging of prostate cancer patients, who are currently staged as pathologic stage 2.

With this method, for the first time, we can directly measure one or more irregularly shaped tumor area(s) on one microscopic slide. Due to limited research conditions, we did not directly incorporate our program into an automatic digital imaging system of the digitally scanned slides. We believe the accuracy will be further improved by directly dotting around the border of cancer areas on digital images. However, based on the same principle, with just a few simple steps, our method can be easily used to acquire objective and accurate tumor volume on a scanned slide or an H&E picture. Furthermore, this method can be also applied to any other solid tumor, whenever an accurate tumor volume/percentage are needed, in particular if multiple separate small tumor foci exist in the same area. Because of its simplicity and ease of application, this mathematical method is an ideal tool to deal with this challenging situation.

Another potential application of our method is for the reconstruction of a histological 3-D tumor image. The current 3-D radiology imaging technique can nicely show the 3-D picture of a prostate mass. However, imaging cannot replace histology. The 3-D mass on radiology imaging does not necessarily represent the size of an actual tumor that we are concerned about.

Some parts of radiologically defined mass may contain benign prostatic hyperplasia (BPH) or be entirely benign. In other words, any suspected imaging has to be confirmed by histology examination, to evaluate the real volume or size of the tumor. Using the currently presented method, we can histologically confirm the cancer borders on serial sections of a prostate gland, with the accurate percentage information on all sections. This allows us to construct a whole 3-D picture of the distribution of the pathologically confirmed tumor in the prostate gland. In addition to Gleason score, grading group, histological type, and molecular information, the accurate data regarding the real size of cancer, its location, distances from margins and adjacent structures should provide clinically significant information for optimal management of individual prostate cancer patients. Future studies comparing the radiologic and pathologic 3-D images of the prostate cancers, may lead to development of management strategies without traumatic surgery.

## Funding

The work was not supported by any commercial fund.

## Availability of Data and Materials

The corresponding author (jchen3@northwell.edu) on reasonable request can provide raw data of results and JavaScript Program.

## Authors' Contributions

John C. conceived the mathematical methods and made the JavaScript program; Jianli C. proposed the clinical question, John, Jianli, MW, and LO wrote the manuscript. All authors read and approved the final manuscript.

## Ethics Approval and Consent to Participate

Not applicable.

## Competing Interests

None.

## References

- <https://www.cancer.org/cancer/prostate-cancer/about/key-statistics.html>
- <https://documents.cap.org/protocols/cp-malegenital-prostate-radicalprostatectomy-19-4041.pdf>
- Morash C, Tey R, Agbassi C, Klotz L, McGowan T, et al. (2015) Active surveillance for the management of localized prostate cancer: guideline recommendations. *Can Urol J Assoc* 9: 171-178.
- Chen RC, Rumble RB, Loblaw DA, Finelli A, Ehdai B, et al. (2016) Active surveillance for the management of localized prostate cancer (Cancer Care Ontario Guideline): American Society of Clinical Oncology practice guideline endorsement. *J Clin Oncol* 34: 2182-2190.
- Sauter G, Steurer S, Clauditz TS, Krech T, Wittmer C, et

- al. (2016) Clinical utility of quantitative gleason grading in prostate biopsies and prostatectomy specimens. *Eur Urol* 69: 592-598.
6. Choy B, Pearce SM, Anderson BB, Shalhav AL, Zagaja G, et al. (2016) Prognostic significance of percentages and architectural types of contemporary Gleason pattern 4 prostate cancer in radical prostatectomy. *Am J Surg Pathol* 40: 1400-1406.
  7. Bismar TA, Lewis JS JR., Vollmer RT, Humphrey PA (2003) Multiple measures of carcinoma extent versus perineural invasion in prostate needle biopsy tissue in prediction of pathologic stage in a screening population. *Am J Surg Pathol* 27: 432-440.
  8. Amin M, Boccon-Gibod L, Egevad L, Epstein JI, Humphrey PA, et al. (2005) Prognostic and predictive factors and reporting of prostate carcinoma in prostate needle biopsy specimens. *Scand J Urol Nephrol* 39: 20-33.
  9. Paner GP, Magi-Galluzzi C, Amin MB, Srigley JR (2014) Adenocarcinoma of the prostate. In: Amin MB, Grignon DJ, Srigley JR, Eble JN, *Urological Pathology*. Lippincott William & Wilkins, Philadelphia, PA, 559-673.
  10. van der Kwast T, Amin M, Billis A, Epstein J, Griffiths D, et al. (2011) The ISUP prostate consensus group. International Society of Urological Pathology (ISUP) Consensus Conference on handling of radical prostatectomy specimens: Working group 2: T2 substaging and prostate cancer volume. *Mod Pathol* 24: 16-25.
  11. Renshaw AA, Richie JP, Loughlin KR, Jiroutek M, Chung A, et al. (1999) Maximum diameter of prostatic carcinoma is a simple, inexpensive, and independent predictor of prostate-specific antigen failure in radical prostatectomy specimens: Validation in a cohort of 434 patients. *Am J Clin Pathol* 111: 641-644.
  12. Perera M, Lawrentschuk N, Bolton D, Clouston D (2014) Comparison of contemporary methods for estimating prostate tumour volume in pathological specimens. *BJU Int* 113: 29-34.

## Supplemental Data:

JavaScript program steps based on Gauss's area calculation formula

```
<!DOCTYPE html>
<html lang = "en">
<head>
<meta charset="utf-8">
<title> FileUpload </title>
<style>
#hello{
  position: relative
}
</style>
</head>
<body>
<input type = "file" />
<button type = "button" />
  <!--determines connection to JS -->
<!-- Button class -->
<p id="showCount"></p>
<p id="showCount1"></p>
<p id="showCount2"></p>
<div class="btn-group">
<button onclick= "CountFun()">Result</button>
<button onclick= "CountFun1()">Tumor</button>
<button onclick= "CountFun2()">Tissue</button>
<button onclick= "LastConnect()">EndPath</button>
</div>

<canvas id="myCanvas" width="1200" height="800" style="border:1px solid #000000;"></canvas>
<script src="script.js"> charset = "utf-8" </script>

<script>
var shapes = 10;
  function LastConnect(){
    shapes+=10;
var a = pointsarray.length;
var canvas = document.getElementById("myCanvas");
var context = canvas.getContext('2d');
context.beginPath();
  context.moveTo(pointsarray[a-1].x, pointsarray[a-1].y);
```

```
context.lineTo(pointsarray[0].x, pointsarray[0].y);
context.stroke();

var i; var A = 0; var B = 0;
pointsarray.push(pointsarray[0]);
for(i=0; i<=pointsarray.length-2; i++)
{
  let a1 = pointsarray[i].x;
  let a2 = pointsarray[i+1].x;
  let b1 = pointsarray[i].y;
  let b2 = pointsarray[i+1].y;
  A+= a1*b2;
  B+= a2*b1;
}
if(boolean)
{
  var element = document.createElement("P");
  element.innerHTML = "Tumor " +cnt1+ " Area: " +0.5*Math.abs(B-A);
  element.id = "hello";
  element.style.left = "500px";
  element.style.top = parseInt("shapes", 10)+"px";
  document.body.appendChild(element);
  totaltumor += 0.5*Math.abs(B-A);
}
if(!boolean)
{
  var element = document.createElement("P");
  element.innerHTML = "Tissue " +cnt2+ " Area: " + 0.5*Math.abs(B-A);
  element.id = "hello";
  element.style.left = "500px";
  element.style.top = parseInt("shapes", 10)+"px";
  document.body.appendChild(element);
  totaltissue += 0.5*Math.abs(B-A);
}

pointsarray.splice(0, pointsarray.length);
}
</script>

</body>
```

```
</html>
```

Javascript:

```
//uploads image. Stretching will not change ratio.
const input = document.querySelector("input[type='file']")
const img = new Image ();
input.addEventListener('change', function(){
  const reader = new FileReader()
  reader.onload = function() {

    img.src = reader.result;
    var c = document.getElementById("myCanvas");
    var ctx = c.getContext('2d');
    ctx.drawImage(img, 0, 0, 1200, 800);
  }
  reader.readAsDataURL(input.files[0])
}, false)

//Button Functions
var cnt=0;
function CountFun(){
  cnt++
  if(cnt>=1)
  {
    window.alert("Tumor/Tissue Percentage: " + (100*totaltumor)/(totaltissue) + "%");
  }
}

var cnt1=0; var boolean;
function CountFun1(){
  cnt1++
  var divData=document.getElementById("showCount1");
  divData.innerHTML="Tumor Pieces: (+cnt1 +)";//this part has been edited
  boolean = true;
}

var cnt2=0;
function CountFun2(){
```

```
cnt2++
var divData=document.getElementById("showCount2");
divData.innerHTML="Tissue Pieces: (" +cnt2 +")";//this part has been edited
boolean = false;
}

var totaltumor = 0;
var totaltissue = 0;

//Button Funtions
var pointsarray = [];

function addpoint(event) {
var rect = myCanvas.getBoundingClientRect();
var docx = (event.clientX - rect.left);
var docy = (event.clientY - rect.top);
var canvas = document.getElementById("myCanvas");
var context = canvas.getContext('2d');
context.beginPath();
context.arc(docx, docy, 5, 0, 2*Math.PI);
context.stroke();
context.fillStyle = "#89F401";
context.fill();

var point = {x: docx, y: docy};
pointsarray.push(point);

if(pointsarray.length > 1)
{
let a = pointsarray.length;
context.beginPath();
context.moveTo(pointsarray[a-2].x, pointsarray[a-2].y);
context.lineTo(pointsarray[a-1].x, pointsarray[a-1].y);
context.stroke();
}
}
myCanvas.addEventListener('click', addpoint);
```